



# AV-Ace II

Pulse Oximeter

# AV-Ace II

## Pulse Oximeter

### Application

Pulse Oximeter can detect SpO<sub>2</sub> and pulse rate through finger. The AV-Ace II handheld pulse oximetry is suitable for hospital's operation room, clinic section office, out-patient department, sickroom, emergency treatment, and the recovery & health care organizations, or in the family nursing and in the process of transporting patients. Etc

- Portable and Light weight
- Easy and simple operation
- Low power consumption
- Display SpO<sub>2</sub>, PR, Bargraph and Plethysmogram
- Low battery indicator
- 2.8 inch TFT LCD
- Adjustable Lightness
- Friendly menu
- Suitable for Adult,
- Pediatric, Neonate



### SpO<sub>2</sub>

Range : 35%-100%  
Resolution : 1%  
Accuracy : 80%-100% ±2%  
70%-80% ±3%  
<70% unspecified

### PR

Range : 25-250bpm  
Resolution : 1%  
Accuracy : ±2%bpm

### Temperature (optional)

Range : 20-50°C  
Resolution : 1%  
Accuracy : ±0.2°C



### AVEUS INSTRUMENTS Sp. z.o.o.

17/19/133 Domaniewska Street,  
02-672 Warsaw - Poland

ul. Domaniewska 17/19 lok.  
133 02-672 Warszawa - Poland

info@aveusinstruments.com.pl  
www.aveusinstruments.com.pl