


2 Years ☆ ☆
Warranty



Light Weight
& Ultra Portable
Design



Built-in
Rechargeable
Battery



Harmless
Infrared
Light




5 level
Adjustable
Brightness



3400 mAh
Built-in
Battery



Reproduction of this catalogue or Pulsemed® logo are liable for legal action. Images shown are for representation purposes only. Engineering Representation in Sizes & Color are approximate. Technical illustrations might be revised for continuous improvements as per the changing industry demands & trends.




VEINUS

Portable Veinfinder



TECHNICAL SPECIFICATIONS

- Detection Distance : 15 - 25 cm
- Accuracy : ± 0.3 mm
- Noise : ≈ 20 dB
- Battery : Rechargeable Li-Ion
- Battery Capacity : 3,400 mAh
- Size : 20x6.2xh 5.5 cm
- Weight : 280 g
- Working time per charge : 175 Mins
- Light : harmless infrared light
- Display: 5 levels adjustable brightness
- Power supply :
- Output : 5V 2.0A DC
- Input : 100V-240V AC @ 50Hz-60Hz



PACKING LIST

- Aluminium Case
- Surgical skin pen
- Tourniquet
- Calibration card
- Charging cable
- Adapter
- User manual
- Tabletop Mount Stand
- Warranty Card
- Charging Adapter
- Main Unit Protect Cover

